

## Phusion DNA polymerase

### Cat. No:

FPLF004.0250

FPLF004.0500

### Contents:

Component	100RXN	200RXN
Phusion DNA poly. 5 U/μl	50 μl	100 μl
MgCl <sub>2</sub> Solution 25 mM	500 μl	1 ml
5X Phusion Buffer MgCl <sub>2</sub> free	500 μl	1 ml

### Description:

Phusion DNA Polymerase is a chimeric Pfu which has a DNA binding protein at the N-terminal portion of the gene. This enzyme keeps significant activity after exposure to 99 °C or repeated exposure to 98 °C with more processivity and extension rate than Pfu DNA polymerase. It catalyzes the polymerization of nucleotides into duplex DNA in the 5'→3' direction, resulting in blunt-ended PCR products **without 3'-dA overhangs**.

Phusion DNA Polymerase exhibits 3'→5' exonuclease (proofreading) activity that enables the polymerase to correct the mis-incorporation of nucleotide, and lacks 5'→3' exonuclease activity. It is suitable for PCR and primer extension reaction that requires high fidelity when the PCR fragment is relatively **higher than 3 kb**. The enzyme exhibits 3'→5' proofreading activity, resulting in over 20-fold higher PCR fidelity than possible with Taq DNA Polymerases.

### Kit storage:

This kit should be stored at -20 °C. Under this condition reagents are stable for two years from the date of production.

Component	Volume	Final Conc.
5X Reaction Buffer	4 μL	1X
MgCl <sub>2</sub> Solution 25 mM	1.2 μL	1.5 mM
40 mM dNTPs Mix (10 mM each)	0.4 μl	0.2 mM
Upstream Primer (10 pmol/ μL)	1 μl	0.5 pmoles/μl
Downstream Primer (10 pmol/ μL)	1 μl	0.5 pmoles/μl
Template DNA	Variable	10 fg~1 μg
PCR grade water	Variable	-
Phusion DNA poly. 5 U/μl	0.25 μl	-
Total Volume	20μl	-

### Protocol:

- 1) Thaw 5X reaction buffer, dNTP mixture.
- 2) Mix the master mix thoroughly and dispense appropriate volumes into PCR tubes or plates.
- 3) Add templates DNA to the individual PCR tubes or wells containing the master mix.
- 4) Program the PCR machine according to the program outlined.

Cycle	Time	Temp °C
1	4 min	95
30-35	30 sec	94
	30 sec	57
	60 sec	72
1	5 min	72

### Note:

- \*Longer extension time makes nonspecific bands
- \*Extension rate for this enzyme is near 3000 bp/min.

### Agarose gel Electrophoresis:

Run the total 5-7 μL of PCR products alongside 3μL DNA marker on a 2% agarose gel containing Green viewer DNA safe stain.