

## Pfu DNA polymerase

### Cat. No:

FPLF002.0250

FPLF002.0500

### Contents:

Components	100 RXN	200 RXN
Pfu DNA poly. 5 U/μl	250 U	500 U
MgCl <sub>2</sub> Solution 25 mM	0.5 ml	1 ml
5X Pfu Buffer MgCl <sub>2</sub> free	0.5 ml	1 ml

### Kit storage:

This kit should be stored at -20 °C. Under this condition reagents are stable for two years from the date of production.

### Description:

Pfu DNA Polymerase is a thermostable enzyme with a molecular weight of 90 kDa. It catalyzes the polymerization of nucleotides into duplex DNA in the 5'→3' direction, resulting in blunt-ended PCR products **without 3'-dA overhangs**. Pfu DNA Polymerase exhibits 3'→5' exonuclease (proofreading) activity that enables the polymerase to correct the mis-incorporation of nucleotide, and lacks 5'→3' exonuclease activity. It is suitable for PCR and primer extension reaction that requires high fidelity when the PCR fragment is relatively **shorter than 3 kbp**. The enzyme exhibits 3'→5' proofreading activity, resulting in over 10-fold higher PCR fidelity than possible with Taq DNA Polymerases.

### Protocol:

- 1) Thaw 5X reaction buffer, dNTP mixture.
- 2) Mix the master mix thoroughly and dispense appropriate volumes into PCR tubes or plates.
- 3) Add templates DNA to the individual PCR tubes or wells containing the master mix.

Component	Volume	Final Conc.
5X Reaction Buffer	4 μl	1X
MgCl <sub>2</sub> Solution 25 mM	1.6 μl	2 mM
40 mM dNTPs Mix (10 mM each)	0.4 μl	0.2 mM
Upstream Primer (10 pmol/ μL)	1 μl	0.5 pmoles/μl
Downstream Primer (10 pmol/ μL)	1 μl	0.5 pmoles/μl
Template DNA	Variable	10 fg~1 μg
PCR grade water	Variable	-
Pfu DNA poly. (5 units/μl)	0.25 μl	-
<b>Total Volume</b>	<b>20μl</b>	<b>-</b>

- 4) Program the PCR machine according to the program outlined.

Cycle	Time	Temp °C
1	4 min	95

30-35	30 sec	94
	30 sec	57
	60 sec	72
1	5 min	72

### Note:

- \* Longer extension time makes nonspecific bands
- \* Extension rate for this enzyme is near 500bp/min.

### Agarose gel Electrophoresis:

Run the total 5-7 μL of PCR products alongside 3 μL DNA marker on a 2% agarose gel containing DNA safe stain.

### Disclaimers and Addresses:

This product is for **Research use only** and should only be used by trained professionals.