

M-MLV Reverse Transcriptase

Cat. No: FPLF011.0100

Store at -20 °C

Contents:

Component	Volume
M-MLV Reverse Transcriptase (100U/ μl)	0.1 ml
5X RT Buffer	0.5 ml

Description:

This a genetically modified RNA-dependent DNA polymerase requiring a DNA primer and an RNA template to synthesize a complementary DNA strand. Thermo-resistant H Minus M-MuLV Reverse Transcriptase has no RNase H activity. Therefore, degradation of RNA does not occur during first strand cDNA synthesis, resulting in higher yields of full-length cDNA from long templates compared to other reverse transcriptases. Thermo-resistant H Minus M-MuLV Reverse Transcriptase maintains activity over a wide temperature range (42-52°C) which makes it an ideal tool for reverse transcription of RNAs having a high degree of secondary structure.

Kit storage:

This kit should be stored at -20 °C. Under this condition reagents are stable for one years from the date of production.

Protocol:

- 1) Mix the template RNA (total RNA or Poly (A) mRNA) and the primer in RNase-free tube as below table. Optimal reaction conditions, such as amount of RNA and primers, may vary and must be individually determined. Random hexamer or oligo (dT) 16 or specific primers could be used as primer.

* If you use RNase inhibitor

Concentration of template RNA and primer

Template RNA	Total RNA or Poly(A)+ mRNA	10 ng~5 μg
	Oligo (dT)16	5 ng~0.5 μg
Primer	or Random hexamer	1-2 μL
		1 μL
DEPC-treated water		Up to 12 μL (11 μL*)

- 2) Incubate the mixture at 65 °C for 5 min and chill on crash ice and add the reagent as follow:

Component	Volume (μL)
5X RT Buffer	4
RNase Inhibitor 20 U/ul (optional)	1
10 mM dNTP Mix	2
Thermo-Resistant RT	2

- 3) Mix by pipetting gently up and down (total reaction volume 20 μL).

- 4) Incubate 10 min at 25 °C (omit this for Oligo dt).
- 5) Incubate 60 min at 47 °C.
- 6) Stop the reaction by heating at 70 °C for 10 minutes. Chill on ice.

Disclaimers and Addresses:

This product is for **Research use only** and should only be used by trained professionals.