

Add-Probe RT-PCR Master (2x Conc.)

Research Use Only

Product Code

74201

Component

1. Add-Probe RT-PCR Master (2x conc.) 1.0 ml

Storage Condition

Store at -20°C

Description

Add-Probe RT-PCR Master (2x) provides sensitive and easy-to-use components which contain all the reagents for first strand cDNA synthesis and PCR reaction in one-tube by using TaqMan®, and is designed for high sensitivity and specificity on various real-time instruments.

Quality Control.

The performance of Add-Probe RT-PCR Master (2x) is tested in an RT and PCR one-tube reaction using human total RNA with specific primers and TaqMan® probes. The sensitivity of the kit is verified by the detection of GAPDH and Actin transcript in 10 pg total RNA after 30 cycles.

Storage and Stability

Add-Probe RT-PCR Master is stable for 2 years when stored in a constant temperature freezer at less than -20°C

Made in KOREA

This product was manufactured through ISO 9001 & 13485 system.

Nucleic Acid Amplification Protocol

1. Add the following components to a thin-walled PCR tube:

| | |
|------------------------------------|-------------|
| Nuclease-Free H ₂ O | x µl |
| Add-Probe RT-PCR Master (2x conc.) | 10.0 µl |
| Forward primer (10 µM) | 0.25~2.0 µl |
| Reverse primer (10 µM) | 0.25~2.0 µl |
| TaqMan Probe (10 µM) | 0.25~2.0 µl |
| (Optional) 50x ROX dye | x µl |
| RNA template | x µl |
| Total reaction volume | 20 µl |

* Recommendation for template RNA concentration in a 20 µl reaction volume

- 1) total RNA: 100 fg ~ 1 µg
- 2) mRNA: 10 fg ~ 1 µg

2. PCR cycling

| | |
|---------------------------------|-----------------------------------|
| cDNA synthesis | 50°C, 20 min |
| Initial denaturation | 95°C, 10 min |
| PCR cycling (30 – 40 cycles) | 95°C, 10 sec 60°C, 30 - 60 sec |

[Note] 50x ROX dye

ROX dye can be included in the reaction to normalize the fluorescent reporter signal, for instruments which are compatible with that option.

50x ROX is a 25 µM concentration. Use the following table to determine the amount of ROX to use with a particular instrument.

| Instrument | Final ROX concentration |
|---|-------------------------|
| AB 7000, 7300, 7700, 7900HT, 7900 Fast, StepOne and StepOnePlus | 500 nM |
| AB 7500, 7500 Fast, Stratagene Mx3000P, Mx3005P and Mx4000 | 50 nM |