



KIAGENE FANAVAR

KIASMID PLUS

Plasmid DNA Isolation Kit (solution based)

Cat. No:

Eco Kit FPKT020.0025

Mini Kit FPKT020.0050

Midi Kit FPKT020.0100

Maxi Kit FPKT020.0200

Contents:

Components	25 preps	50 preps	100 preps	200 preps
PSB1 Buffer (blue)	2.5 ml	5 ml	10 ml	20 ml
PSB2 Buffer (orange)	5 ml	10 ml	20 ml	40 ml
PSB3 Buffer (yellow)	3.75 ml	7.5 ml	15 ml	30 ml
PSB4 Buffer (green)	12.5 ml	25 ml	50 ml	100 ml
PEB	1 ml	2 ml	4 ml	8 ml
Rnase A 10mg/ml	100ul	200ul	400ul	800ul

Kit storage:

⚠ This kit should be stored at room temperature.

⚠ RNase A should be stored at -20 °C

If properly stored, all kit components are stable until the expiration date printed on the label.

⚠ After adding RNaseA to PSB1 buffer Store the PSB1 at +2 to +8°C

Additional Equipment and Reagent required

- Absolute ethanol
- Ethanol 70%
- Standard tabletop microcentrifuge capable of 13,000 x g centrifugal force
- Microcentrifuge tubes, 1.5 ml, sterile

Application:

Isolation of up to 15µg purified plasmid DNA from bacterial cultures, which may be used directly in downstream applications such as restriction enzyme digestion, PCR, cloning, sequencing, in vitro transcription, or labeling reactions. The procedure is optimized to achieve reliable results within 90 min.

Handling Requirements and Safety Information

⚠ All solutions are clear and should not be used when precipitates have formed. Warm the solutions at +15 to +25°C or in a 37°C water bath until the precipitates have dissolved.

⚠ PSB4 contain Phenol which is an irritant and toxic.

⚠ Do not allow PSB4 Buffers to touch your skin, eyes, or mucous membranes. If contact does occur, wash the affected area immediately with large amounts of water.

⚠ If you spill the reagent, dilute the spill with water before wiping it up.

⚠ Do not use any modified ethanol.

⚠ Do not pool reagents from different lot numbers.

⚠ Immediately after usage, close all bottles in order to avoid leakage, varying buffer concentrations or buffer conditions.

⚠ After first opening store all bottles in an upright position.

⚠ Do not allow the PSB4 Buffers to mix with sodium hypochlorite found in commercial bleach solutions. This mixture can produce a highly toxic gas.

⚠ Wear protective disposable gloves, laboratory coats and eye protection, when handling samples and kit reagents.

⚠ Do not contaminate the reagents with bacteria, virus, or nucleases. Use disposable pipets and nuclease free pipet tips only, to remove aliquots from reagent bottles.

preparation procedure:

Working Solution Preparation

⚠ add RNase A to PSB1 Buffer completely Dissolve the RNase A by vortexing. store the PSB1 buffer at 4 °C.

Protocol

1. Sample preparation

Pellet the bacterial cells from **1-5 ml of E. coli culture** (OD600= 1.5-5 per ml) by centrifuge in **8000g for 1 min**. Discard the supernatant (Do not use more highly concentrated sample)

2. Add 100 µl chilled PSB1 Buffer (with RNase A) to the bacterial pellet. Resuspend the bacterial pellet by **vortex for 30 s**.

3. Add 200 µl PSB2 Buffer to the mixture. Mix gently by **inverting the tube 5 to 10 times (Do not vortex)**. **Incubated at room temperature for 1 min** (Do not incubate more).

4. Add 150 µl chilled PSB3 Buffer. Mix gently by **inverting the tube 5 to 10 times (Do not vortex)**. **Incubated on ice for 5 min** (The solution should become cloudy and a flocculent precipitate should form).

5. Centrifuge at Maximum speed for 5 min at +4°C and transfer entire supernatant into a new 1.5 ml microtube.

Contact and Support: If you have questions or experience problems with Kiagene Fanavar products, please contact our Technical Support staff. Our scientists are committed to providing rapid and effective help. Website: www.kiagene.ir Email: Techsupport@kiagene.ir Tel: 02191010809

6. Add 500 µl PSB4 and Centrifuge at **Maximum speed for 2 min at +4°C**

7. Transfer upper solution layer into a new microtube.
Add 2-Fold Absolute ethanol and incubate 2 min at RT.

8. Centrifuge at Maximum speed for 5 min at +4°C. (a thin and white pellet should be seen in bottom of microtube) Discard supernatant carefully.

9. Add 1ml Ethanol 70% and invert several. Centrifuge at **Maximum speed for 2 min at +4°C.**

10. Discard supernatant carefully. to remove residual ethanol, **incubate microtube at +55°C on hot plate for 5-10 min.**

11. Add 50 µl PEB Buffer or sterile DDW on pellet and pipetage gently.

12. Check **5-10 µl of extracted plasmid DNA** by 1% agarose gel electrophoresis. Store Plasmid DNA at -20°C